

# Release Notes

---

## V 4.1.1 (release on 2026-06-23)

### Bug fixes

- Fixed an issue where ARCSolver exits with a license error while a valid pressure serial number license is present
- Fixed an issue where ARCSolver considers only one serial number, if more than one licenses with serial numbers are present

## Overview of previous releases of ARCSolver

| Date       | Version | Changes   |
|------------|---------|---|
| 2026-04-30 | 4.1.0   | <p><b>Features</b></p> <ul style="list-style-type: none"> <li>• <b>Link the ARCSolver license to the serial number of the measurement hardware</b></li> <li>• <b>Use the integrated license tool to create and activate a license for another offline ARCSolver installation</b></li> </ul> <p><b>Improvements/Changes</b></p> <ul style="list-style-type: none"> <li>• Added support for Linux on ARM64 architectures</li> <li>• Added integration example code for Python 3</li> <li>• More flexible timestamps for ultrasound diameter analysis</li> <li>• Updated SOUPS: <ul style="list-style-type: none"> <li>• OpenSSL to 3.5.5</li> </ul> </li> </ul> <p><b>Bug fixes</b></p> <ul style="list-style-type: none"> <li>• Fixed an issue with the rpath of libcrypto.3.dylib on macOS (x86_64 only)</li> </ul> |

| Date       | Version | Changes  |
|------------|---------|--|
| 2025-09-30 | 4.0.0   | <p><b>Features</b></p> <ul style="list-style-type: none"> <li>• <b>Process instantaneous arterial diameter curves recorded with ultrasound</b> (new input modality)</li> </ul> <p><b>Improvements/Changes</b></p> <ul style="list-style-type: none"> <li>• Activate the software license without internet connection</li> <li>• Control feature usage with different licenses</li> <li>• New interface for encrypted signal transfer</li> <li>• Numerical optimization to improve inter-platform compatibility</li> <li>• Updated SOUPS: <ul style="list-style-type: none"> <li>• Qt to 6.8.3</li> <li>• OpenSSL to 3.5.0</li> <li>• SoftwareKey to 5.25.1.0 (Desktop)</li> <li>• SoftwareKey to 5.23.4.1 (Android)</li> </ul> </li> </ul> <p><b>Bug fixes</b></p> <ul style="list-style-type: none"> <li>• Do not write files into program folder</li> <li>• Improved description of input signals band-pass filter</li> <li>• Update integration manual to better cover error cases</li> </ul> |
| 2024-02-29 | 3.0.0   | <p><b>Features</b></p> <ul style="list-style-type: none"> <li>• <b>first release of ARCSolver under MDR</b></li> <li>• ARCSolver is now stand-alone microservice</li> </ul> <p><b>Improvements</b></p> <ul style="list-style-type: none"> <li>• Added license control for post market surveillance</li> <li>• Removed deprecated input parameter "sigtype"</li> <li>• Separated input parameter "flag" into "calibration" &amp; "huckDetection" for better readability</li> </ul> <p><b>Bug fixes</b></p> <ul style="list-style-type: none"> <li>• change to 64-bit data types to improve inter-platform compatibility</li> <li>• Minor bugfix in calculation of AIX based on 3rd derivative</li> <li>• Abort when more than 30 pulse waves are detected (to prevent memory overflow)</li> <li>• Abort if hr &lt; 40 or hr &gt; 140 during calculation to stay within defined physiological boundaries</li> </ul>  |
| 2019-03-19 | 2.0.2   | <ul style="list-style-type: none"> <li>• Not for general use!</li> <li>• SV calculation optimization</li> </ul>  |
| 2017-07-06 | 2.0.1   | <ul style="list-style-type: none"> <li>• check for poor end-diastolic signal quality ("Huck", probably due to wrong recording pressure) set error code to 1 when poor quality detected</li> <li>• enable / disable end-diastolic signal quality detection ("Huck") using map_flag</li> </ul>   |
| 2015-12-09 | 2.0.0   | <ul style="list-style-type: none"> <li>• <b>Removed debug output from developer version for production version 1.8.a</b></li> </ul>  |

| Date       | Version   | Changes  |
|------------|-----------|--|
| 2015-12-03 | 1.8.a     | <ul style="list-style-type: none"> <li>• DEVELOPER VERSION WITH SENSITIVE DEBUG ISSUE!</li> <li>• Pulse wave analysis identical to 1.7.2</li> <li>• Conversion to encrypted input data</li> <li>• Instead of the parameter signal of the type double[1000], which represents ADC values, raw data of the type double[1533 .. 2033] are now passed, which represent encrypted and possibly packed byte values.</li> </ul> |
| 2015-09-08 | 1.7.2.S.1 | <ul style="list-style-type: none"> <li>• <b>Release for Schiller AG, Delphi based.</b></li> </ul>  |
| 2014-04-22 | 1.7.2     | <ul style="list-style-type: none"> <li>• <b>Bugfix: In exceptional cases, the check as to whether the Alx is too high was not carried out</b></li> <li>• <b>Bugfix: Improving the central curve at the inflection point sometimes caused the calculation to crash</b></li> <li>• <b>Bug fix for evaluation with mean pressure calibration</b></li> </ul>   |
| 2013-12-18 | 1.7.0     | <ul style="list-style-type: none"> <li>• Improvements for central pulse wave form estimation</li> </ul>  |
| 2012-02-22 | 1.6.5     | <ul style="list-style-type: none"> <li>• Improvement of the Alx calculation for large Alx values</li> <li>• Improvement of optimization from rel. 1.6.4</li> </ul>   |
| 2012-02-10 | 1.6.4     | <ul style="list-style-type: none"> <li>• Correction for flattening curves (e.g. nitro administration)</li> <li>• Removal of calculation for individual curves</li> </ul>   |
| 2012-01-25 | 1.6.3     | <ul style="list-style-type: none"> <li>• <b>Optimization of the calculation of the stroke volume when specifying the mean pressure</b></li> <li>• <b>Correction in the PWV calculation for children</b></li> </ul>   |
| 2011-10-04 | 1.6.2     | <ul style="list-style-type: none"> <li>• Improved Vascular Age calculation</li> </ul>  |
| 2011-07-15 | 1.6.1     | <ul style="list-style-type: none"> <li>• Bugfix: Optional calibration was not applied to the radial pulse curve. Was corrected.</li> </ul>   |
| 2011-04-08 | 1.6.0     | <ul style="list-style-type: none"> <li>• Optional transfer and calibration of the central pulse curve with the peripheral mean pressure</li> </ul>   |
| 2011-04-01 | 1.5.2     | <ul style="list-style-type: none"> <li>• <b>Optimization in the PWV and Alx calculation</b></li> </ul>   |
| 2011-02-17 | 1.5.1     | <ul style="list-style-type: none"> <li>• Improved Alx calculation and plausibility check</li> </ul>  |
| 2010-09-21 | 1.5.0     | <ul style="list-style-type: none"> <li>• New Feature: Vascular Age calculation</li> </ul>  |
| 2010-08-18 | 1.4.5     | <ul style="list-style-type: none"> <li>• First release integrated into commercial 24h blood pressure monitor Mobil-O-Graph (IEM)</li> </ul>  |